IFU-OTH-0004

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Revisão / Revision: 0
Data / Date: 07/09/2023

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PRODUCT DESCRIPTION

Flexy NiTi orthodontic archwires are made from resilient Nitinol shape memory alloys, which are heat-treated to deliver constant force loads that meet clinical expectations. In addition, the archwires feature variations in composition with the addition of chemical elements that provide additions of elastic, malleable and aesthetic characteristics, providing the patient with comfort, quality, safety and reliability.

When used in conjunction with brackets, tubes and accessories fixed to the tooth surface, the archwires aim to assist in movement or force generation, helping to establish the mechanical-biological system, and thus enable the application of force to the teeth, which results in induced tooth movement, correcting malocclusion.

The Flexy NiTi archwires have differentiating features such as the layout which provides a more functional arch form for occlusion, greater arch extension which makes it possible to include second molars, and also the use of stops, marking of the midline, differentiating into three lines for the upper arch and one line for the lower arch, making it easier to identify them.

The Flexy line of archwires is made up of alloys and formats that allow them to be used at different stages of treatment.

- Flexy NiTi Super Elastic - TrueForm and ALX formats;
- Flexy NiTi Super Elastic, Esthetic Rhodium - TrueForm and ALX formats;
- Flexy NiTi Thermal 35 °C - TrueForm and ALX formats
- Flexy NiTi Thermal 35 °C, Esthetic Rhodium - TrueForm and ALX formats
- Flexy NiTi Copper - TrueForm and ALX formats;
- Flexy NiTi Copper, Esthetic Rhodium - TrueForm and ALX formats;
- Flexy NiTi Reverse Curve - TrueForm and ALX formats;
- Flexy Blue-Ti - TrueForm and ALX formats.

The upper and lower archwires are available in the following sizes: Ø.012"; Ø.014"; Ø.016"; Ø.018"; Ø.020"; .014"x.025"; .016"x.016"; .016"x.022"; .017"x.025"; .018"x.025"; .019"x.025"; .021"x.025.

Flexy NiTi Super Elastic archwires have been developed for use in the stages of orthodontic treatment involving alignment, rotation and leveling. Flexy NiTi Thermal 35 °C archwires are thermally activated and have a lower load-deflection than conventional NiTi archwires. Activations below room temperature and body temperature facilitate the archwires' adaptation in the oral cavity.

The Flexy NiTi Reverse Curve archwire group is characterized by promoting the extrusion of posterior teeth and buccalization of incisors, thus helping to correct deep overbite.

The inclusion of 6% copper resulted in Flexy NiTi Copper, which has better shape memory properties and greater resistance to permanent deformation when compared to conventional alloys, guaranteeing excellent control of tooth movements.

The aesthetic requirement is achieved by the Esthetic Rhodium line, which is found in the Flexy NiTi Super Elastic and Flexy NiTi Copper lines and is characterized by rhodium plating that promotes an opacity in the arches and a whitish tone, making it closer to tooth colouring, guaranteeing greater aesthetics.

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
Finally, the Flexy Blue-Ti line has a smoother and more polished surface than conventional NiTi alloys, which reduces friction and makes it possible to use them in mechanical space closure by sliding, while preserving their superelastic characteristics, allowing them to be used in the initial stages of orthodontic treatment, alignment and leveling.

Flexy NiTi archwires should be chosen according to the specific characteristics of the clinical case and are identified according to the drawing below:

INTENDED USE

Flexy NiTi archwires fall within the range of Nickel and Titanium archwires and have the mechanical properties that meet the needs of a working archwires, which are: low modulus of elasticity, high flexibility, minimal plastic deformation, resistance to fatigue, good resistance to corrosion, production of light and constant forces.

These excellent mechanical properties combined with NiTi's shape memory and super elasticity make Flexy NiTi archwires ideal for the initial stages of orthodontic treatment. During these stages, dental elements show irregularities in alignment and leveling, to a greater or lesser degree of severity, and Flexy NiTi archwires can undergo deflections to adjust to the brackets on the teeth, releasing light and constant forces, resulting in bone response and physiological tooth movements.

The shape memory and super elastic properties are responsible for the light and constant forces released by the deflections of the archwires, even when unloaded, reducing the number of wire changes and providing less discomfort for patients.

NiTi Flexy archwires and their variations are used together with brackets, tubes and orthodontic bands in the initial stages of aligning and leveling Class I, Class II and Class III sagittal malocclusions, preparing the dental arches for the subsequent stages of orthodontic mechanics.

In transverse malocclusions, characterized by posterior crossbites, Flexy NiTi archwires can be used after the orthopedic maxillary expansion phase, to promote alignment and leveling of the dental elements, enabling the use in sequential phases of stiffer and heavier archwires, such as stainless-steel alloys.

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In the case of vertical problems involving deep overbite and open bite, Flexy NiTi archwires can enable alignment and leveling of the upper and lower incisors before intrusive and extrusive mechanics are carried out. This allows the problems to be corrected with lighter and more constant forces, ensuring physiological movements with less demand on the anchorage units.

Inter-arch and intra-arch problems such as mild and moderate crowding, tooth rotation and inadequate axial inclination can be corrected with the use of Flexy NiTi archwires, but associated with ways of opening spaces, such as interproximal wear or tooth expansion and protrusion mechanics.

PRINCIPLES OF FUNCTIONALITY

Flexy NiTi archwires together with brackets, tubes and bands make up the mechanical-biological system and make it possible to apply the light and constant forces necessary for induced tooth movement, in a controlled and physiological manner and with an adequate bone response.

The movements that can be achieved are:

- 1- Vertical intrusion and extrusion.
- 2- Rotation; and
- 3- Mesiodistal angulations.

INSTRUCTIONS FOR USE

For the arches to show their full potential, we recommend the following technical rigor in their use:

Application / Installation of the arches:

1. The professional must select the appropriate archwire according to the needs of the treatment;
2. Measure the length of the archwire to be used in the oral cavity beforehand on a study model;
3. Cut off the excess archwire;
4. Prophylaxis and polishing of the brackets and, above all, cleaning of the accessory channels;
5. Selecting the archwire section consistent with the alignment and leveling irregularities of the dental elements, in order to respect the mechanical properties of the alloy;
6. Guide the archwire into the oral cavity with Weingart pliers or needle-nose instruments and insert it into the mesial entrance of the molar tube;
7. Fit the archwire into the bracket channels and start tying the archwire to the bracket fins;
8. Fix the archwire with elastic or metal ligatures on conventional brackets or fix the archwire by closing the clips on self-ligating brackets;
9. Tie the archwire to the flaps of the bracket of the tooth positioned mesially to the tooth with the most suitable positioning, and then make a larger deformation towards the bracket of this element and tie it to the adjacent teeth;

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IFU-OTH-0004

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Revisão / Revision: **0**
Data / Date: 07/09/2023

10. Cut off the excess archwire that appears distal to the tube, using a pair of distal cutters which have the ability to cut and fix the wire segment.
11. If you wish to make a distal bend, first break off the end of the archwire before initially inserting the archwire into the bracket channels, and then bend it using Weingart pliers or pliers suitable for distal bending.
12. Check the ligatures.

Wire Removal:

1. Remove the elastic or metal ligatures in conventional brackets or open the clips in the case of self-ligating brackets;
2. Remove the archwire by traction in a mesial direction.

STORAGE CONDITIONS

Keep away from sunlight.

Keep dry.

Do not use if packaging is damaged.

CONTRAINDICATIONS

Orthometric only supplies its products to qualified professionals. It is the orthodontist's sole responsibility to identify conditions that may be contraindicated for treatment, such as:

- Patients with poor oral hygiene;
- Patients unable to cooperate with treatment;
- Allergy to any of the orthodontic appliance materials;
- Diseases or other pre-existing conditions that may hinder orthodontic treatment;
- Existing bone resorption;
- Existing decalcification of the tooth enamel.

SIDE EFFECTS

It is the responsibility of the qualified professional to identify any possibility of risk of harm and/or contraindications that may arise during treatment, to inform the patient of the possibility of any unwanted side effects.

During treatment, unwanted side effects may include: periodontal complications, allergic reactions, difficulties in maintaining oral hygiene, discomfort and pain.

WARNINGS

The product is supplied in non-sterile condition, however, it has proven biocompatibility and no cleaning, disinfection or sterilization measures are required, provided the packaging has not been tampered with. If the product gets into contact with surfaces or substances other than those intended for use, it is recommended that it be discarded.

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Heat treatment is not recommended.

In the case of aesthetic archwires, to avoid damaging the aesthetic layer, replace the intraoral archwire within 30 days.

Be careful when handling the archwire with pliers or tweezers and install accessories such as stops and hooks with caution.

Heating the archwire above 200°C can induce fracture and, in the case of aesthetic archwires, deteriorate the aesthetic layer.

Excessive stress and sharp bends can fracture the archwires.

Frequent exposure to fluoride-containing substances can cause the archwire to become brittle.

For distal cutting in the oral environment, ensure that the pliers retain the tip of the archwire.

This product should not be used in an MRI environment, as it may cause interference in the images during the examination.

This product contains Nickel. A small percentage of the population is allergic to this metal. If allergic reactions occur, advise the patient to seek medical advice.

This product is not reusable, as reuse can cause cross-infection and loss of mechanical properties due to natural wear. The product must not be reused or reprocessed.

Do not use the product if its packaging is damaged or open, in order to minimize possible dangers related to the transmission of microorganisms or possible product counterfeiting.

This product should only be used by qualified orthodontic professionals. Use by people without the necessary technical knowledge could lead to unwanted tooth movement, bone loss, tooth loss or even irreversible damage. It is up to the qualified professional to advise the patient on proper oral hygiene, in order to avoid the appearance of plaque or diseases such as gingivitis, periodontitis or even endocarditis.

Considering that the product is used in contact with mucous membranes and body fluids, it is recommended that the professional applies good clinical practice for the disposal of medical products in accordance with current health regulations.


SIMBOLOGY

	Do not re-use	<table border="1"> <tr> <td data-bbox="853 1758 933 1825">EC</td> <td data-bbox="933 1758 1013 1825">REP</td> </tr> </table>	EC	REP	Authorized representative in the European Community
EC	REP				

IFU-OTH-0004

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Revisão / Revision: 0
Data / Date: 07/09/2023

	Consult instructions for use		European conformity
	Medical device		Caution
	Catalogue number		Non-sterile
	Batch code		Contains Nickel
	Manufacturer		Date of manufacture
	Use-by date		Keep away from sunlight
	Keep dry		Do not use if package is damaged

ATTENTION: Orthometric cannot be held responsible for the results obtained by applying the product without the necessary precautions or by failing to observe the warnings and/or by application by professionals who are not qualified to use the products.