IFU-OTH-0005

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Revisão / Revision: 2
Data / Date: 18/12/2025

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IFU-OTH-0005

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Revisão / Revision: 2
Data / Date: 18/12/2025

INSTRUCTIONS FOR USE
English version

Technical name: Molar Distalizer

Product name: CLASS II CORRECTOR - ICERAM CLEAR DISTALIZER

ANVISA nº: 80328819003, 80328819011

Important: To consult the Instructions for Use, access the digital file on the website: www.orthometric.com.br/instrucoes_de_uso or via the QR code provided on the packaging.

To obtain a printed copy of the Instructions for Use, please contact our Customer Service (SAC) at phone number 0800 770 4045, mobile +55 (14) 99690-1718 or e-mail sac@orthometric.com.br.

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Product Description

ICeram aesthetic distalizers are orthodontic devices designed to promote the distal movement of posterior teeth – primarily in the maxilla – with the goal of improving the anteroposterior relationship between the dental arches. They act by distalizing the posterior teeth as a unit, while maintaining patient comfort and an enhanced esthetic appearance.

Available in seven different sizes: 14.0 mm, 16.0 mm, 18.0 mm, 20.0 mm, 23.0 mm, 25.0 mm, and 27.0 mm, the ICeram aesthetic distalizers are manufactured from engineering-grade materials with esthetic properties, making them a suitable option for aligner-based treatments. Their main features include: ease of installation; universal design suitable for use on either the right or left side of the arch; esthetic and comfortable.

This is a fixed appliance used in the treatment of symmetrical or asymmetrical Class II malocclusions, intended for the distalization of maxillary molars, premolars, and canines. It is well accepted by patients due to its simple use and discreet design, eliminating the need for attachments on all teeth. Additionally, it enables independent, straight-line distalization of the posterior segment without resistance from the curved anterior segment, which accelerates anteroposterior correction and eliminates unwanted deflection forces.

Intended Use

The ICeram Distalizer is designed for use in the treatment of Class II malocclusions, enabling correction through distalization of the posterior segment by promoting distal tipping and rotation of the molars. To apply the required force, it must be used in combination with Class II elastics. The device may be used with fixed appliances in the mandibular arch or with aligners. The elastic attachment point should be positioned on the mesial portion of the accessory and on a metal or ceramic button, or a tube bonded to the buccal surface of the mandibular first molar.

When used with fixed appliances in the lower arch, reinforcement of anchorage with a lingual arch is recommended. When used in combination with aligners in the mandibular arch, the aligner must be trimmed in the area where the esthetic button is bonded, ensuring proper adaptation of the aligner tray.

Indications for Use

The ICeram Distalizer is part of a biomechanical system used in conjunction with brackets, tubes, buttons, and orthodontic archwires or aligners, enabling the application of the force required to achieve controlled tooth movement.

The movements that can be achieved include distalization of the maxillary posterior teeth and distalization of the mandibular posterior teeth.

Instructions for Use

1. Size Selection

To determine the appropriate distalizer size, follow the steps below:

- i. Use the measuring ruler provided with the product.
- ii. When used for movements beginning at the canines, measure the distance between the mesial third of the maxillary first molar and the center of the buccal surface of the maxillary canine.
- iii. When used for movements beginning at the second premolar, measure the distance between the mesial third of the maxillary first molar and the center of the mesial portion of the maxillary second premolar.
- iv. If the patient's measurement does not correspond exactly to any of the available product sizes, select the next larger option (e.g., if the patient's measurement is 17.0 mm, use the 18.0 mm distalizer).
- v. Measurements must be taken on both the right and left sides of the patient to identify possible asymmetry requirements.

2. Installation of the ICeram Distalizer

Once the appropriate size has been selected, proceed with the bonding protocol, following the steps below:

- i. Isolate the working field using plastic cheek retractors and relative isolation with cotton rolls, keeping a high-power saliva ejector in place to prevent contamination of the bonding areas with saliva.
- ii. Identify and delimit the specific areas intended for accessory bonding.
- iii. Perform prophylaxis with pumice, water, and a Robinson brush to remove contaminants.
- iv. Rinse the cleaned areas thoroughly with water.
- v. Etch only the enamel surface where the accessory will be bonded using 37% phosphoric acid gel for 15 to 30 seconds. The purpose of acid etching is to create micro-porosity for mechanical interlocking of the bonding resin.
- vi. Rinse the etched areas thoroughly, keeping the high-power suction in place to avoid contamination of the surfaces.
- vii. Dry the etched enamel surfaces with an air jet free of water or oil contamination.
- viii. Apply primer to the enamel surface in the bonding area using a brush without applying pressure.
- ix. Light-cure the primer.
- x. The application of a Plastic Primer, Universal Adhesive, and/or Bond Enhancer (Assure or Assure Plus, Reliance) on the esthetic base to be bonded is recommended.
- xi. Apply a layer of resin to both the esthetic and metal bases.
- xii. Position the accessory on the enamel surface using a bonding tweezer, applying light pressure.
- xiii. Adapt the excess resin around the base of the accessory, creating additional reinforcement for adhesion.
- xiv. Light-cure for 3 seconds on each base.
- xv. Insert the Class II elastic, size ¼" medium or 3/16". It is recommended to begin with 170 g of force during the first month and increase to 230 g until correction of the molar relationship is achieved.

3. Removal of the Distalizer

To remove the Distalizer, use a straight or angled bracket-removal plier, according to professional preference and adaptation, applying traction and shear force.

The technical removal procedure involves the following steps:

- i. Remove any excess resin that may be surrounding the accessory base. However, if bonding was performed with proper technique, minimal resin should be present around the orthodontic accessory.
- ii. Use Straight or Angled Bracket Removal Pliers, or Weingart pliers, applying traction and shear force to deform the accessory and release it from the enamel surface through failure at the accessory–resin interface. The choice of instrument depends on the clinician's familiarity and preference.
- iii. Remove adhesive using scraping pliers with carbide-tipped or titanium tips. Excess resin may also be removed with a bracket-removal plier.
- iv. Remove remaining composite using a tapered multilaminated bur made of tungsten carbide, used in a low-speed contra-angle handpiece, with 9 blades. A high-speed 24-blade bur may also be used. These burs feature reverse blades to preserve the dental enamel.
- v. After removing resin remnants, finish the enamel surface using a Shofu Ceramisté polishing rubber. Final polishing is performed with rubber cups and Spanish white in water, or with cloth wheels in a contra-angle handpiece and diamond paste.

Composition

Hinge – AISI 630 (17-4 ph), 30%

Arm – Polyamide type 12, 70%

Storage Conditions

Store in a clean and dry place, protected from intense light, at room temperature, and away from oxidizing agents or contaminants. Do not place objects and/or weights on top of the product.

The manufacturer recommends the single use of this product. It should be disposed of correctly in accordance with the information contained in the Manual of Dental Services Risk Prevention and Control, Technology in Health Services, 2006 (available on the ANVISA website).

Warnings and Symbolology

- It is essential that only a qualified specialist installs the product on the patient. In addition to performing the clinical examination, the professional must have the necessary knowledge to systematically describe the patient's dentofacial characteristics, list priority problems, and outline possible solutions for each issue, considering potential favorable or unfavorable interactions for each proposed approach.
- It is the orthodontist's responsibility to determine the maximum expected benefit from orthodontic therapy and the minimum acceptable loss, or to calculate the benefit-risk ratio for each patient, including psychosomatic factors, to ensure proper case management.
- Patients must be instructed regarding oral hygiene, maintenance, and dietary care. The following should be avoided: ice, coconut, hard and/or sticky candies, caramel, nougat, peanut brittle, corn on the cob, popcorn, and sugarcane. Biting hard objects such as pens and pencils causes similar issues.
- Advise the patient NOT TO CHEW OR BITE HARD FOODS during orthodontic treatment to prevent premature debonding.
- This product must not be used in a Magnetic Resonance Imaging (MRI) environment, as it may cause image interference during the exam. It is the orthodontist's responsibility to inform the patient about this contraindication and arrange removal of the appliance if needed.
- The high concentration of bacteria and fungi present in the oral cavity can accelerate the corrosion of metal orthodontic components. The pH of the oral environment also greatly influences the corrosion process; therefore, a diet rich in acidic foods – such as excessive consumption of acidic soft drinks – should be avoided as much as possible during treatment.
- Do not touch the components with ultrasonic scaler tips; this may cause fractures or debonding.
- The product is provided non-sterile; however, it has proven biocompatibility, and no cleaning, disinfection, or sterilization is required as long as the packaging has not been opened. If the product comes into contact with surfaces or substances outside its intended use prior to application, it must be discarded.
- The product must not be reused or reprocessed.
- Do not use the product if its packaging is damaged or opened, in order to minimize risks related to microorganism transmission or potential product counterfeiting.
- It is the responsibility of the qualified professional to inform patients of the risk of aspiration or swallowing

in the event of breakage/detachment of components. The clinician must not underestimate potential immediate or delayed complications arising from such incidents and must properly inform the patient and/or legal guardians about inherent risks and the appropriate actions to be taken.

- This product is intended exclusively for use by qualified orthodontic professionals. Use by individuals without the necessary technical knowledge may result in undesirable tooth movement, bone loss, tooth loss, fenestration, gingival recession, or root dehiscence.
- It is the responsibility of the qualified professional to instruct the patient on proper oral hygiene to prevent bacterial plaque, calculus, gingival inflammation (gingivitis), periodontitis, or even endocarditis. As the product is used in contact with mucosa and bodily fluids, the professional must follow clinical best practices for disposal of medical devices in accordance with applicable health regulations.

ATTENTION: Orthometric cannot be held responsible for the results obtained by applying the product without the necessary precautions or failing to observe the warnings.

Contraindications

Orthometric provides its products exclusively to qualified professionals. It is the orthodontist's full responsibility to identify conditions that may contraindicate treatment, such as:

- Patients with poor oral hygiene.
- Patients unable to cooperate with treatment.
- Patients with an oral environment prone to enamel demineralization.

	Do not re-use		Non-sterile
	Consult instructions for use		Contains Nickel
	Medical device		Caution
	Catalogue number		Date of manufacture
	Batch code		Keep away from sunlight
	Manufacturer		Do not use if package is damaged

IFU-OTH-0005

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Revisão / Revision: **2**
Data / Date: 18/12/2025

	Use-by date		Keep dry
---	-------------	---	----------