IFU-OTH-0001

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- Telemóvel: +351 967 417 654

Product description

The metal bracket is made of stainless steel and has an anatomical design that provides a good fit to the dental elements, providing the patient with comfort, quality, safety, and reliability.

The base has pins that help retain the bracket to the adhesive used for bonding for better adhesion of the adhesive to the bracket.

The base design provides marginal sealing, preventing resin leakage.

The brackets have gemined t-wing, the groove is .022" x.028" and have chamfered edges in order to reduce friction.

When fixed to the buccal surface of the tooth, the brackets aim to house the intraoral orthodontic arches and interact with other orthodontic accessories, helping to establish the mechanical-biological system, thus enabling the application of force to the teeth, which results in induced tooth movement, correcting malocclusion. In some cases, brackets can serve as supports for the use of cantilevers, which increase their versatility in tooth movements.

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Intended use

Brackets are used when orthodontic mechanics are required to correct the various malocclusions, sagittal, vertical and transverse, inter-arch and intra-arch, taking into account the appropriate diagnosis in all cases. The following malocclusions can be treated using brackets:

I- Sagittal problems:

- 1) Angle Class I malocclusion;
- 2) Angle Class II malocclusion;
- 3) Angle Class III malocclusion;
- 4) Mild mandibular prognathism (in compensatory treatment) and in severe cases can be used in pre-surgical orthodontic preparation and in the post-surgical phase to provide adequate occlusion;
- 5) Mild mandibular retrognathism (in compensatory treatment in adults, and in patients in the active growth phase, making up the mandibular propulsion system) and in severe cases can be used in pre-surgical orthodontic preparation and in the post-surgical phase to provide occlusion adequacy;
- 6) Mild maxillary protrusion (compensatory treatment in adults, and in growing patients after the orthopedic phase) and in severe cases can be used in pre-surgical orthodontic preparation and in the post-surgical phase to provide occlusion adequacy;
- 7) Mild maxillary retrusion (compensatory treatment in adults, and in growing patients after the orthopedic phase) and in severe cases can be used in pre-surgical orthodontic preparation and in the post-surgical phase to provide occlusion adjustment;
- 8) Anterior crossbite (by correcting the malocclusion through buccal inclination of the upper incisors and/or lingual inclination of the lower incisors);

II- Transverse problems:

- 1) Posterior crossbite (in the post-orthopedic expansion phase or as a compensatory treatment for the buccal inclination of postero-superior teeth) and also after orthognathic surgery has been performed and the post-surgical containment phase has elapsed.

III- Vertical problems:

- 1) Deep overbite;
- 2) Anterior open bite (in orthodontic or orthodontic-surgical treatment) and in severe cases they can be used in pre-surgical orthodontic preparation and in the post-surgical phase to adjust the occlusion.

IV- Inter-arch and intra-arch problems:

- 1) Dental crowding;
- 2) Midline deviations;
- 3) Dental positioning asymmetries.

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Revisão / Revision: **3**
Data / Date: 27/05/2024

How to use:


For proper bonding of brackets, the following sequence of procedures must be adopted and are summarized as:

1. prophylaxis;
2. acid etching of the enamel;
3. application of the primer followed by bonding the bracket to the tooth.

The instructions for bonding are described below:

DIRECT BONDING

1. Isolation of the work field with plastic cheek retractors and relative isolation with cotton rolls, keeping the high-power saliva sucker, to avoid contamination of the gluing areas with the sage; Delimitation of specific areas for gluing accessories;
2. Perform prophylaxis with a pumice stone, water, and a Robson brush to remove contaminants;
3. Wash the sanitized areas with plenty of water;
4. Promote the conditioning of the enamel surface, only the area that will receive the bracket, with the use of 37% phosphorus acid gel applied for 15 to 30 seconds. This acid etching procedure aims to create microporosity for the mechanical imbrication of the bonding resin;
5. Wash the conditioned areas intensively, keeping the high power suction in position, to avoid contamination of these surfaces.
6. Drying of conditioned enamel surfaces with a jet of air free from contamination with water or oil;

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

Revisão / Revision: 3
Data / Date: 27/05/2024

7. Application of fluid resin or primers on the enamel surface referring to the bonding area, with a brush and without exerting pressure;

8. Photopolymerization;
9. Keep the bracket fixed between the ends of the bonding forceps;
10. Application of a thin layer of primers throughout the bracket base;
11. Positioning the bracket on the enamel surface with the use of collage forceps, applying light pressure;
12. Conference of occlusal retraction or X;
13. Removal of excess resin;
14. Light cure for 3 seconds per bracket.

II- Indirect bonding

- 1) Selection of the brackets of interest to the professional;
- 2) Choosing the trays;

- 3) Prophylaxis of the upper and lower dental arches to eliminate the presence of any contaminating physical agent from the enamel surfaces;
- 4) Molding the upper and lower dental arches with irreversible hydrocolloid (alginate);
- 5) Obtaining type IV stone plaster models or digital models;
- 6) Delimitation of the long axes of the dental elements buccally and lingually, and demarcation of the EV points and X bonding heights on the plaster models;
- 7) Applying plaster insulation to the plaster surfaces in the places that will receive the resin brackets;
- 8) Leave to dry for 30 seconds;
- 9) Bonding the attachments with gold quality orthodontic resin, removing excess resin and checking the positioning of the attachments;
- 10) Light curing: follow the resin manufacturer's instructions;
- 11) Obtaining the transfer trays.

Direct and indirect bonding techniques may vary depending on the professional who performs them.

Due to the large number of steps involved in the bracket bonding process, it is very difficult to control the moisture absorbed in the process, which can lead to contamination of the enamel surface and, consequently, failure to adhere to the tooth. The presence of moisture and/or contamination during the bonding process is the biggest cause of bracket failure, delaying and increasing treatment costs.

Brackets can undergo changes to their surface structure when exposed to the adverse conditions present in the oral environment. In the oral cavity, metal biodegradation usually occurs through electrochemical disruption. Although electrolytic attack seems to be the major cause of disruption, bacteria and their cellular degradation products, and selective interactions with gases such as oxygen and carbon dioxide can contribute to the degradation of materials in the oral cavity. Organic acids and enzymes in particular can also affect the various metals and the pH of the environment where the devices are installed, exerting a major influence on the corrosive process.

The medical product in question may only be fitted to a patient by a dental surgeon specializing in orthodontics.

Bracket removal

Once the treatment objectives have been achieved and the finishing touches have been made, the brackets are removed with the adhesives without damaging the tooth enamel structure. It is recommended that suitable, high-quality instruments are used to avoid discomfort for the patient and the risk of injury to soft tissue and dental structures.

The technical removal procedures involve the following steps:

- 1) Removal of excess resin that may be surrounding the base of the brackets, however, if the bonding has been carried out with technical care, the minimum possible amount of resinous material around the orthodontic accessory is observed;
- 2) Use of Straight Bracket Remover Pliers or Weingart pliers, both of which are used to perform mesiodistal compression of the bracket fins, which causes deformation of the accessory and its release from the enamel surface by breaking the bracket-resin interface. The instruments used at this stage will depend on the

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Revisão / Revision: 3
Data / Date: 27/05/2024

professional's familiarity with the products.

- 3) Remove the adhesive using scraper pliers with Widia or titanium tips. Excess resin can also be removed with a pair of pliers;
- 4) Removal of composite residues with a multilaminate conical trunk drill made of tungsten carbide, adapted to a 9-blade low rotation contra-angle. You can also use the high-rotation 24-blade drill, code CF 375R. These drills have inverted blades to preserve tooth enamel; and
- 5) After removing the resin residue, the surface of the enamel should be finished with a Shofu Ceramisté polishing rubber. The final polishing is carried out with rubber cups and Spain white in water or cloth wheels in a contra-angle and diamond paste.

Composition:

Metal Bracket - 17-4 PH steel (Designation UNS S17400)

Storage Conditions:

In intact and closed packaging.

This product does not require temperature, humidity, or light control during storage. The recommendation is to keep it in a clean, dry place.

The manufacturer recommends the single use of this product. It should be disposed of correctly in accordance with the information contained in the Manual of Dental Services Risk Prevention and Control, Technology in Health Services, 2006 (available on the ANVISA website).

Warnings and Symbology:

- It is essential that only the specialist installs the product on the patient, since in addition to the clinical examination, the professional must have the necessary knowledge to make a systematic description of the patient's dentofacial features, list the priority problems, describe possible solutions for each of the problems, considering the possible favorable or unfavorable interactions for each of the proposed attempts.
- It is also the orthodontist's duty to determine the maximum expected gain from orthodontic therapy and the minimum loss, or to calculate the benefit versus risk ratio for each patient, with both risk and benefit including psychosomatic factors to determine the best management for each patient.
- Patients should be instructed in oral hygiene, maintenance, and food care. The following should be avoided: ice, coconut, hard and/or sticky sweets, caramels, nougats, jellybeans, corn, popcorn and sugar cane. Biting hard objects such as pens and pencils has the same effect. Chewing gum and toffee often loosen the bands or rings that have been cemented in place.
- Advise patients NOT to chew or bite on hard foods during orthodontic treatment to prevent premature detachment. This product should not be used in an MRI environment, as it may cause interference in the images during the examination. It is up to the orthodontic professional to warn the patient of this contraindication and arrange for the appliance to be removed if necessary.
- The high concentration of bacteria and fungi present in the mouth can accelerate the corrosion of

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Revisão / Revision: **3**
Data / Date: 27/05/2024

metal orthodontic appliances. The pH of the oral environment also has a major influence on the corrosive process, so a diet rich in acidic foods, such as excessive intake of soft drinks with acidic pHs, should be avoided as much as possible during treatment.

- Do not touch brackets with ultrasonic tips: this can lead to fractures and detachment.
- The product is supplied in non-sterile condition, however it has proven biocompatibility and no cleaning, disinfection or sterilization measures are required, provided the packaging has not been tampered with. If the product comes into contact with surfaces or substances other than those intended for use, it is recommended that it be discarded.
- The product must not be reused or reprocessed. Do not use the product if its packaging is damaged or open, to minimize possible dangers related to the transmission of microorganisms or possible product counterfeiting.
- It is the responsibility of qualified professionals to advise their patients of the risk of aspiration and swallowing in the event of breakage/ detachment of parts, not to underestimate the possible immediate or delayed complications resulting from this type of accident and to adequately inform the patient and/or their legal guardians of the risks inherent in the event and the measures to be taken.
- This product should only be used by qualified orthodontic professionals.
- Use by people without the necessary technical knowledge could lead to unwanted tooth movement, bone loss, tooth loss, fenestration, gingival recession, or root dehiscence.
- It is the responsibility of the qualified professional to advise the patient on correct oral hygiene, in order to avoid the appearance of plaque and tartar or diseases such as gingival inflammation (gingivitis), periodontitis or even endocarditis. Considering that the product is used in contact with mucous membranes and body fluids, it is recommended that professionals apply good clinical practices for the disposal of medical products in accordance with current health regulations.

ATTENTION: Orthometric cannot be held responsible for the results obtained by applying the product without the necessary precautions or failing to observe the warnings.

Contraindications


Orthometric only supplies its products to qualified professionals. It is the orthodontist's sole responsibility to identify conditions that may be contraindicated for treatment, such as:

- Patients with poor oral hygiene.
- Patients unable to cooperate with treatment.
- Patients with an oral environment conducive to the demineralization of tooth enamel.

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Revisão / Revision: **3**
Data / Date: 27/05/2024

	Do not re-use		Authorized representative in the European Community
	Consult instructions for use		European conformity
	Medical device		Caution
	Catalogue number		Non-sterile
	Batch code		Contains Nickel
	Manufacturer		Date of manufacture
	Use-by date		Keep away from sunlight
	Keep dry		Do not use if package is damaged