IFU-OTH-0011

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Revisão / Revision: 0
Data / Date: 15/01/2024

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Revisão / Revision: 0
Data / Date: 15/01/2024

**INSTRUCTIONS FOR
USE**

Technical Name: Components for Orthodontics

Product Name: ORTHOMETRIC TUBE PREMIUM SERIES

ANVISA n°: 80328810065

Technical Name: Components for Orthodontics

Product Name: ORTHOMETRIC TUBE ADVANCED SERIES

ANVISA n°: 80328810065

Technical Name: Metallic Brackets

Product Name: ORTHOMETRIC TUBE CAPELOZZA SYSTEM

ANVISA n°: 80328810074

Technical Name: Metallic Brackets

Product Name: ORTHOMETRIC TUBE M JANSON

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The tubes are accessories made of stainless steel and they are welded welded to the molar bands or glued directly to molar teeth, or occasionally to premolars. Depending on its design, it can be called simple, when it only has one fitting for the orthodontic arch; double tube, which is characterized by fitting the internal arch of the headgear and the main arch; triple, when it has three fittings, one for the orthodontic arch, another for auxiliary devices and the third for the internal arch of the headgear; double grooved tubes, which serve as fittings for rectangular arches; and lingual or palatal tube for adapting the transpalatal bar.

The Orthometric Premium Series, Advanced Series, Capelozza System and Mjanson tube are made by metal injection molding (MIM) which makes possible to bring together the requirements of quality, technology and

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The Orthometric Premium Series tubes are manufactured according to the Roth prescription, with a 0.022" x 0.028" groove, being: 1) simple non-convertible tube with hook for welding and gluing; 2) simple convertible tube with hooks for welding and gluing; 3) round/rectangular convertible double tube for welding; 4) Rectangular convertible double tube for welding and gluing; and 5) round and double rectangular convertible triple tube for welding. For the MBT prescription, the tubes available are: 1) simple non-convertible tube for welding and gluing; 2) convertible single tube for welding and gluing; and 3) round and rectangular convertible double tube for welding

The Advanced Series, Capelozza and MJanson tubes feature an anatomical design that establishes a good adaptation to the coronal surfaces of the dental elements of interest, providing the patient with comfort, quality, safety and reliability.

Advanced Series tubes are available with a 0.022" x 0.028" groove in simple non-convertible models, with hooks, for gluing or welding in Roth and MBT prescriptions; and simple non-convertible model for gluing with or without hooks in Edgewise and Ricketts prescriptions

Capelozza tubes are available with a 0.022" x 0.028" groove in the simple convertible, simple non-convertible, and double convertible models for gluing in the Capelozza I prescription; and simple convertible for gluing in the Capelozza II prescription.

MJanson tubes are available with a 0.022" x 0.028" groove in simple convertible, simple non-convertible and simple mini non-convertible second molar models for bonding in the MBT prescription.

The Advanced series, Capelozza and MJanson tubes are characterized by having: 1) superior notch that allows better positioning of accessory; 2) permanent engraving U and L, differentiating them respectively for the upper and lower arches; 3) tapered entry to facilitate arch insertion; 4) ease of manipulation with the tweezers with thin, angled tips, which helps in grasping the tube by the tweezers during gluing; 5) anatomical base curved in the mesiodistal and occlusogingival direction, allowing better adaptation to the tooth surface; 6) retention pins on the base ensuring greater adhesion efficiency; and 7) rounded hooks, which ensures greater comfort for the patient, and allows the use of elastics.

The series of Orthometric Premium series, Advanced series, Capelozza and MJanson tubes guarantee the professional a versatility of options that meet clinical needs.

When welded to the band or fixed by means of gluing to the buccal surface of the tooth, the tubes establish the final support of the fixed appliance, and aim to house the intraoral orthodontic arches, helping to establish, together with the brackets, the mechanical-biological system, and in this way, they enable the application of force to the teeth, which results in induced tooth movement, correcting the malocclusion.

In some cases, the double and triple tubes can serve as supports for the use of cantilevers, and when necessary, they can serve as support for the internal arch of the headgear, contributing to the application of force that redirects the growth of the jaw. All applications for tubes increase their versatility in tooth movements.

Intended use

Tubes are used when orthodontic mechanics are required to correct the various malocclusions, sagittal, vertical and transverse, inter-arch and intra-arch, taking into account the appropriate diagnosis in all cases. The following malocclusions can be treated:

I- Sagittal problems:

- 1) Angle Class I malocclusion;
- 2) Angle Class II malocclusion;
- 3) Angle Class III malocclusion;
- 4) Mild mandibular prognathism (in compensatory treatment) and in severe cases can be used in pre-surgical orthodontic preparation and in the post-surgical phase to provide adequate occlusion;
- 5) Mild mandibular retrognathism (in compensatory treatment in adults, and in patients in the active growth phase, making up the mandibular propulsion system) and in severe cases can be used in pre-surgical orthodontic preparation and in the post-surgical phase to provide occlusion adequacy;
- 6) Mild maxillary protrusion (compensatory treatment in adults, and in growing patients after the orthopedic phase) and in severe cases can be used in pre-surgical orthodontic preparation and in the post-surgical phase to provide occlusion adequacy;
- 7) Mild maxillary retrusion (compensatory treatment in adults, and in growing patients after the orthopedic phase) and in severe cases can be used in pre-surgical orthodontic preparation and in the post-surgical phase to provide occlusion adjustment;
- 8) Anterior crossbite (by correcting the malocclusion through buccal inclination of the upper incisors and/or lingual inclination of the lower incisors);

II- Transverse problems:

- 1) Posterior crossbite (in the post-orthopedic expansion phase or as a compensatory treatment for the buccal inclination of postero-superior teeth) and also after orthognathic surgery has been performed and the post-surgical containment phase has elapsed.

III- Vertical problems:

- 1) Deep overbite;
- 2) Anterior open bite (in orthodontic or orthodontic-surgical treatment) and in severe cases they can be used in pre-surgical orthodontic preparation and in the post-surgical phase to adjust the occlusion.

IV- Inter-arch and intra-arch problems:

- 1) Dental crowding;
- 2) Midline deviations;
- 3) Dental positioning asymmetries.

How to use:

The direct and indirect bonding techniques are described below:

DIRECT BONDING



1. Isolation of the work field with plastic cheek retractors and relative isolation with cotton rolls, keeping the high-power saliva sucker, to avoid contamination of the gluing areas with the sage;
2. Delimitation of specific areas for gluing accessories;
3. Perform prophylaxis with a pumice stone, water, and a Robson brush to remove contaminants;
4. Wash the sanitized areas with plenty of water;
5. Promote the conditioning of the enamel surface, only the area that will receive the bracket, with the use of 37% phosphorus acid gel applied for 15 to 30 seconds. This acid etching procedure aims to create microporosity for the mechanical imbrication of the bonding resin;
6. Wash the conditioned areas intensively, keeping the high power suction in position, to avoid contamination of these surfaces.
7. Drying of conditioned enamel surfaces with a jet of air free from contamination with water or oil;
8. Photopolymerization;
9. Keep the bracket fixed between the ends of the bonding forceps;
10. Application of a thin layer of primers throughout the bracket base;
11. Positioning the bracket on the enamel surface with the use of collage forceps, applying light pressure;
12. Conference of occlusal retraction or X;
13. Removal of excess resin;
14. Light cure for 3 seconds per bracket.
15. Application of fluid resin or primers on the enamel surface referring to the bonding area, with a brush and without exerting pressure;

IFU-OTH-0011

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Revisão / Revision: 0
Data / Date: 15/01/2024

II- Indirect bonding

- 1) Selection of the product of interest to the professional;
- 2) Choosing the trays;
- 3) Prophylaxis of the upper and lower dental arches to eliminate the presence of any contaminating physical agent from the enamel surfaces;
- 4) Molding the upper and lower dental arches with irreversible hydrocolloid (alginate);
- 5) Obtaining type IV stone plaster models or digital models;
- 6) Delimitation of the long axes of the dental elements buccally and lingually, and demarcation of the EV points and X bonding heights on the plaster models;
- 7) Applying plaster insulation to the plaster surfaces in the places that will receive the pieces with resin;
- 8) Leave to dry for 30 seconds;
- 9) Bonding the attachments with quality orthodontic resin, removing excess resin and checking the positioning of the attachments;

10) Light curing: follow the resin manufacturer's instructions;

11) Obtaining the transfer trays.

Direct and indirect bonding techniques may vary depending on the professional who performs them.

Due to the large number of steps involved in the bonding process, it is very difficult to control the moisture absorbed in the process, which can lead to contamination of the enamel surface and, consequently, failure to adhere to the tooth. The presence of moisture and/or contamination during the bonding process is the biggest cause of mini tube failure, delaying and increasing treatment costs.

Tubes can undergo changes to their surface structure when exposed to the adverse conditions present in the oral environment. In the oral cavity, metal biodegradation usually occurs through electrochemical disruption. Although electrolytic attack seems to be the major cause of disruption, bacteria and their cellular degradation products, and selective interactions with gases such as oxygen and carbon dioxide can contribute to the degradation of materials in the oral cavity. Organic acids and enzymes can also affect the various metals and the pH of the environment where the devices are installed, exerting a major influence on the corrosive process.

The medical product in question may only be fitted to a patient by a dental surgeon specializing in orthodontics.

Glued Tubes Removal

Once the treatment objectives have been achieved and the finishing touches have been made, the tubes are removed with the adhesives without damaging the tooth enamel structure. It is recommended that suitable, high-quality instruments are used to avoid discomfort for the patient and the risk of injury to soft tissue and dental structures.

The technical removal procedures involve the following steps:

- 1) Removal of excess resin that may be surrounding the base of the tubes, however, if the bonding has been carried out with technical care, the minimum possible amount of resinous material can be observed around the orthodontic accessory;
- 2) Use of Straight Bracket Removal Pliers or Weingart Pliers, both of which are used to remove the accessory and release it from the enamel surface by breaking the part-resin interface. The instruments used at this stage will depend on the professional's familiarity.
- 3) Remove the adhesive using scraper pliers with Wídia or titanium tips. Excess resin can be removed with pliers;
- 4) Removal of composite residues with a multilaminar conical trunk drill made of tungsten carbide, adapted to a counter-angle at low speed with 9 blades. The 24-blade high-speed drill can also be used. The drills have inverted blades to preserve tooth enamel; and
- 5) After removing the resin residue, the surface of the enamel should be finished with a Shofu Ceramisté polishing rubber. The final polishing is carried out with rubber bowls and white Spain in water or cloth wheels in a contra-angle and diamond paste.

Tube Welding

- 1) Selection of tubes of interest;
- 2) Demarcation on the outer buccal surface of the band of a reference line passing through the mesiobuccal cusp of the molar (vertical reference line);

- 3) Put the tube between the electrode tips of the spot welding machine, in such a way that the mesial entry of the tube coincides with the vertical reference line previously defined on the band, and with the tube groove, respecting the X determined for welding;
- 4) Welding at just one point per mesial on the tube fixing base;
- 5) Check the angulation of the tube at the mouth and whether the positioning of the accessory follows the vertical and X reference line;
- 6) Welding of the entire tube, using welding points along the entire length of the fixing base; and
- 7) Cementation of the band with the tube (see step-by-step instructions for banding and cementation).

Removal of he bands

Using band pulling pliers, the removal is carried out, exerting a force to expel the crown band from the dental element, in an occlusal direction. In order for this procedure to be carried out with all criteria and greater comfort for the patient, causing less risk to tooth enamel, the occlusal support platform of the pliers is supported on the occlusal face of the molar, and with the active tip adapted to the cervical region of the pressure band. if in an occlusal direction.

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After removing the band, the cement residues that are adhered to the tooth enamel are eliminated by scraping this material with the wídia tip of the pliers or by using a composite residue removal drill, a multi-bladed trunk drill conical, made of Tungsten Carbide, adapted to a low-speed contra-angle with 9 blades, CB 27, code 75.80.1004.

The final polishing of the tooth surface is carried out with a rubber cup and Spanish white, using a contra-angle and at low speed.

Composition:

Metalic Tubes – 17-4 PH Steel (UNS S17400)

Storage Conditions:

In intact and closed packaging.

This product does not require temperature, humidity, or light control during storage. The recommendation is to keep it in a clean, dry place.

The manufacturer recommends the single use of this product. It should be disposed of correctly in accordance with the information contained in the Manual of Dental Services Risk Prevention and Control, Technology in Health Services, 2006 (available on the ANVISA website).

Warnings and Symbology:

- It is essential that only the specialist installs the product on the patient, since in addition to the clinical examination, the professional must have the necessary knowledge to make a systematic description

IFU-OTH-0011

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Revisão / Revision: **0**
Data / Date: 15/01/2024

of the patient's dentofacial features, list the priority problems, describe possible solutions for each of the problems, considering the possible favorable or unfavorable interactions for each of the proposed attempts.


- It is also the orthodontist's duty to determine the maximum expected gain from orthodontic therapy and the minimum loss, or to calculate the benefit versus risk ratio for each patient, with both risk and benefit including psychosomatic factors to determine the best management for each patient.
- Patients should be instructed in oral hygiene, maintenance, and food care. The following should be avoided: ice, coconut, hard and/or sticky sweets, caramels, nougats, jellybeans, corn, popcorn, and sugar cane. Biting hard objects such as pens and pencils has the same effect. Chewing gum and toffee often loosen the bands or rings that have been cemented in place.
- Advise patients NOT to chew or bite on hard foods during orthodontic treatment to prevent premature detachment. This product should not be used in an MRI environment, as it may cause interference in the images during the examination. It is up to the orthodontic professional to warn the patient of this contraindication and arrange for the appliance to be removed if necessary.
- The high concentration of bacteria and fungi present in the mouth can accelerate the corrosion of metal orthodontic appliances. The pH of the oral environment also has a major influence on the corrosive process, so a diet rich in acidic foods, such as excessive intake of soft drinks with acidic pHs, should be avoided as much as possible during treatment.
- Do not touch the pieces with ultrasonic tips: this can lead to fractures and detachment.
- The product is supplied in non-sterile condition, however it has proven biocompatibility and no cleaning, disinfection or sterilization measures are required, provided the packaging has not been tampered with. If the product comes into contact with surfaces or substances other than those intended for use, it is recommended that it be discarded.
- The product must not be reused or reprocessed. Do not use the product if its packaging is damaged or open, to minimize possible dangers related to the transmission of microorganisms or possible product counterfeiting.
- It is the responsibility of qualified professionals to advise their patients of the risk of aspiration and swallowing in the event of breakage/ detachment of parts, not to underestimate the possible immediate or delayed complications resulting from this type of accident and to adequately inform the patient and/or their legal guardians of the risks inherent in the event and the measures to be taken.
- This product should only be used by qualified orthodontic professionals.
- Use by people without the necessary technical knowledge could lead to unwanted tooth movement, bone loss, tooth loss, fenestration, gingival recession, or root dehiscence.
- It is the responsibility of the qualified professional to advise the patient on correct oral hygiene, in order to avoid the appearance of plaque and tartar or diseases such as gingival inflammation (gingivitis), periodontitis or even endocarditis. Considering that the product is used in contact with mucous membranes and body fluids, it is recommended that professionals apply good clinical practices for the disposal of medical products in accordance with current health regulations.

ATTENTION: Orthometric cannot be held responsible for the results obtained by applying the product without the necessary precautions or failing to observe the warnings.

Contraindications

Orthometric only supplies its products to qualified professionals. It is the orthodontist's sole responsibility to identify conditions that may be contraindicated for treatment, such as:

- Patients with poor oral hygiene.
- Patients unable to cooperate with treatment.
- Patients with an oral environment conducive to the demineralization of tooth enamel.

	Do not re-use		Authorized representative in the European Community
	Consult instructions for use		European conformity
	Medical device		Caution
	Catalogue number		Non-sterile
	Batch code		Contains Nickel
	Manufacturer		Date of manufacture
	Use-by date		Keep away from sunlight
	Keep dry		Do not use if package is damaged