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**Photopolymerizable Adhesive For Orthodontic Bonding**

**EN – Instructions for use**

**Professional use only**

Please, read all the information containing this instruction for use before using the product. Keep it for later consultation until the product is fully consumed, and/or until there is no more interaction between the product and your last patient.

**Product Indication**

Ortholink One Step adhesive, also known as orthodontic composite resin or orthodontic resin cement, was developed for metal, ceramic and polycarbonate orthodontic brackets direct bonding to the dental surface.

**Product Description**

Ortholink – One Step is photopolymerizable and monocomponent. The product is self-adhesive on etched enamel, as the primer, adhesive, and fillers

are combined in their formulation, eliminating, therefore, the adhesive use, providing the professional with clinical ease. Contains the main methacrylic adhesive monomers, such as BisGMA, TEGDMA, and MDP phosphate; being that the polymerization is performed through blue halogen light photoactivation at 400 to 500nm of wavelength or with LED. Barium-Aluminum glass and high-performance Nano Silica distribution can be observed in its composition. Such ingredients in Ortholink - One Step promote a considerable increase in bond strength, adequate working time, allow brackets' easy handling, and avoid their displacement before photopolymerization. The product presents a fluorescent dye that allows its tracking when there is an ultraviolet light incidence, turning orthodontic resin residues identification easier after the brackets installation and/or removal, helping to avoid damage to the enamel.

#### Instruction for use

The instructions for use are described here, however, it is worth mentioning that it is orthodontist's criteria to evaluate the best strategy for the product use depending on each patient's needs for the treatment (the amount of product to be applied, the increments format, and the orthodontic resin in the mouth function time can vary considerably in each treatment). Ortholink One Step is self-adhesive on etched enamel, eliminating, therefore, dental adhesive use in most cases. To improve bonding orthodontic

accessories (Brackets, Tubes, Buttons, etc.) adherence and excellence even more, dental adhesive use is recommended, as well as in cases where the mechanical demand is higher, e.g. overbite, brachyfacial patients, bruxism. 1. Operative field isolation: operative field must be kept free of saliva contamination through an effective isolation. Additionally, for better access to the teeth, a lip retractor use is recommended.

2. Dental preparation: proceed with prophylaxis using oil-free pastes and Robinson brush, or rubber cup. Rinse with plenty of water, and dry the teeth in oil-and moisture-free air.

3. Enamel surface conditioning with 37% of phosphoric acid according to the manufacturer's recommendations. Apply the acid to the bracket bonding area only. Wash with plenty of water, and dry.

4. Dental adhesive may be applied to enamel as per manufacturer's instructions.

NOTE: Ortholink - One Step is a self-adhesive composite, a dental adhesive application is optional and according to each situation needs duly assessed by the orthodontist.

5. Bracket bonding: apply a small amount of Ortholink - One Step directly using a syringe at the bracket base, and gently position it on the previously etched tooth surface. Adjust the bracket in final position, squeezing it firmly to seat it and eliminate excess material. NOTE: The product is ambient light sensitive, in case the final positioning is postponed, it is recommended to work with the reflector turned off momentarily, and/or cover the patient's mouth in order to avoid the material premature polymerization. Clean the applicator tip and replace the syringe cap at the end of your

application.

6. Excess resin removal must be carried out prior to polymerization, excess resin around the orthodontic accessory can be removed with the aid of an explorer probe. Ortholink - One Step, as an auxiliary resource, has a fluorescent (UV) tracker allowing its identification using ultraviolet light sources (between 390 – 410 nm with 250 mW/cm<sup>2</sup> of a minimum power), turning visualization to remove excess resin easier. Such residues removal must be carried out in a careful manner, taking care of not moving the bracket from its correct position and not wearing down and/or damaging the tooth enamel.

NOTE: excesses can be removed, even when an auxiliary UV source prior to polymerization is not used.

7. For metallic and ceramic brackets, use 40 s of polymerization, and 30 s for polycarbonate with photopolymerizer with a minimum power of 600 mW/cm<sup>2</sup>. For metallic and ceramic brackets or polycarbonate with photopolymerizer with a minimum power of 1400 mW/cm<sup>2</sup>, use 20 s of polymerization. A photopolymerizer with a wavelength light between 400 to 500 nm must be used. Perpendicular photopolymerizer to the light always reaching the bracket with the adhesive at a distance of 2 to 3 mm.

8. Archwires can be installed promptly after the brackets bonding.

9. The brackets must be removed using specific instruments after completing the orthodontic treatment. The remaining cement adhered to the tooth surface can be removed using a special multi-laminated drills and fine-grained sandpaper discs, using UV light to identify resinous remnants presence and prevent damage to the enamel. Although the extra fluorescence of the cement makes it easily identified

under UV light, the professional must always be careful when removing product residues. 10. After, the teeth should be polished using felt discs and extra fine-grained polishing paste.

#### Presentation form

Presentation: Pack containing 1 4g syringe of the product and instructions for use for professionals.

#### Basic Composition

Ingredients: Methacrylic monomers (such as BisGMA and TEGDMA), phosphate monomers (such as MDP); photoinitiator, NaF, co-initiator, stabilizers, inorganic Barium-Aluminum fillers, Silica, silane and pigment.

#### Expiration date

Printed on product packaging.

#### Storage conditions

Ortholink - One Step must be stored at temperatures between 5°C to 27°C/ 41°F to 80.6°F. The product must be kept in its original packaging. Keep the product in a cool place, keeping the packaging always well closed. Protect from direct sunlight. Do not store the material next to other eugenol-containing products as this may inhibit resin cement proper polymerization. Do not freeze the product. Protect the product from external contamination. Do not expose the product to intense light or high temperatures.

#### Handling Conditions

The product is intended for professional use only. Read the instructions for use carefully. Do not use the product if it is out of its

expiration date.

For product disposal, please, follow your country's legislation.

After use, clean the instruments with alcohol.

After each use, replace/close the cap.

When handling the product, wear protective gloves.

If the eyes or other oral tissues contact the product by accident, rinse them with large amounts of water.

If necessary, traces of the non-polymerized product can be cleaned with the aid of alcohol.

#### Warnings

Do not use the product if it is out of its expiration date.

Keep out of the reach of children.

#### Precautions

Professional use only

Follow the manufacturer's instructions correctly.

Avoid contact of the non-polymerized product with the skin, mucous membranes and eyes. In case of contact with eyes, rinse them immediately with plenty of water.

Avoid its use in patients with a history of allergy to methacrylate-based composite resins. In case of allergy, discontinue its use and, if necessary, refer the patient to medical assistance.

Suitable protective goggles must be used, when using a photopolymerizer unit.

Light sensitive. Always keep the product closed when it is not in use. The contents must be used quickly after extruding the product out of the syringe.

Care should be taken when applying the product to contact lenses, veneers and prosthetic crowns.

Do not apply the product to thin or

compromised crowns, lenses or veneers. For porcelain pieces, porcelain primer must be used.

Do not store the product close to eugenol-based materials as these may inhibit the cement/adhesive proper curing.

Photopolymerizer equipment must frequently have its power measured with a radiometer.

#### Contraindications and adverse effects

The product contains methacrylate monomers that may cause mild allergic or irritating reactions in these substances-sensitive patients. In case of allergic reactions to the product, discontinue use.

As it is an exclusive use product on enamel, no side effects are expected if the product is used as recommended.

This material was manufactured for dental use only and must be handled according to the instructions for use. The manufacturer is not responsible for damage caused by other uses or incorrect handling. Furthermore, it is mandatory for the user to prove, before using and under his or her responsibility, this material compatibility with the intended use, mainly when these instructions for use do not indicate its type of use. Data descriptions do not constitute any type of guarantee and, therefore, does not present a binding nature.

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